


easy mount[®]

Automatic sample mount and retrieval for the mar dtb





Overview

The **easymount** is an addition to the *marresearch desktop beamline (dtb)*. Its purpose is to safely mount (and retrieve) samples onto (from) the goniometer head. The operation is supported from both a local motor control unit and from the software. The following motors and utilities are added to a standard *dtb*:

(A) Motorized cryo actuator: moves the cryo nozzle out of the way during sample load or unload.

(B) Motorized x,y,z-translations for the goniometer head.

(C) Two additional motors for the sample loading and unloading mechanism.

(D) Illumination unit for well-defined illumination conditions around the sample (not visible).

(E) Extended local motor control unit with buttons for x,y,z and Phi-motions as well as for activating the load/unload mechanism (not shown).

Operation

The UNLOAD operation comprises the following steps:

- (1)** The Chi axis moves up into a vertical position.
- (2)** A vial holding device swings out to an easily accessible position and waits for a vial filled with LN2 to be inserted by the user.
- (3)** The vial holder transfers the vial to the goniometer head.
- (4)** The cap is taken off the head and placed on the vial.
- (5)** The vial holder transfers the vial to a position where it can easily be grabbed by the user.

The LOAD operation reverses the UNLOAD operation:

- (6)** The user inserts a vial with a cap into the vial holder.
- (7)** The vial holder transfers the vial to the goniometer head (blue arrow) and leaves the cap on the goniometer head.
- (8)** The Chi axis moves back down (red arrow).

marresearch

Marresearch GmbH
Hans-Böckler-Ring 17
D-22851 Norderstedt
Germany

Phone: +49 40 529 884- 0
Fax: +49 40 529 884-20
info@marresearch.com
www.marresearch.com

marUSA

marUSA, Inc.
1880 Oak Avenue
Evanston, IL 60201
USA

Toll free: 1-877-MAR-XRAY
From outside U.S.: 1-847-869-1548
Fax: 1-847-869-1587
info@mar-usa.com
www.mar-usa.com