

Comparison of the GeniX ^{3D} Cu HF beam delivery system with the GeniX Cu VHF system

Claudio Klein, Marresearch GmbH, Norderstedt

Introduction



The *mar_{ppX}* system is a very compact yet powerful complete X-ray setup suitable for both protein and small molecule data collection as well as other types of experiments. The *mar_{ppX}* system consists of a *mar₃₄₅* image plate detector, a *mar_{dtb}* goniometer system and the new 2010 model of the GeniX ^{3D} Cu HF generator operated at 30 W (50 kV / 0.6 mA).

In this study we compare the new GeniX ^{3D} Cu HF system with the previous GeniX Cu VHF generator operated at 50 W (50 kV / 1.0 mA). The new system features an entirely new tube setup including a newly designed and highly optimized X-ray optic. Thanks to a number of major improvements the new generator delivers more X-ray photons onto the sample despite of running at lower power (30 W instead of 50W).

In order to assess the quality of the new system, several small molecules of different types and sizes have been collected. Here we present a selection of typical results.



Technical Data

	New GeniX	Old GeniX
Beam at crystal position The picture show an area of 0.5 x 0.5. The color coding reflects the measured intensity of the beam.		
Distance source - optic center	8 cm	14 cm
Distance optic center - focus	30 cm	39 cm
Spot size in focus	175 μm x 175 μm	190 μm x 190 μm
Flux within 100 μm at focus	9 x 10 ⁷ photons/sec	5 x 10 ⁷ photons/sec
Divergence (mrad @ FWHM)	~ 6 x 6	5.4 x 5.4

Data collection and processing

The data for all crystals were collected on the same *mar345 dtb* detector system. When exchanging the generators, the beam was carefully realigned. The crystals needed to be taken off the goniometer but they were remounted in a position very close to where they have been mounted before. Data were collected at:

- 90 mm distance crystal-to-detector giving 1.5 Ang. resolution at the edge
- 5 deg./image
- 72 images (i.e. 360 deg. of data)
- constant exposure time (typically 40s/image)

Data were processed using the *automar* program suite.

	Anil		Bitartrate1		Bitartrate2	
Crystal						
Space group	P 1		P 2 ₁ 2 ₁ 2 ₁		P 2 ₁ 2 ₁ 2 ₁	
Unit cell axes a, b, c	a=9.7 b=13.5 c=16.7		a=7.6 b=7.8 c=11.1		a=7.6 b=7.8 c=11.1	
Mosaicity	0.6°		0.6°		0.3°	
Size of crystal W x H x L	350 x 180 x 500		120 x 120 x 2000		30 x 60 x 500	
Optics	NEW	Old	NEW	Old	NEW	Old
# reflections	4166	4194	1287	1316	1417	1387
Multiplicity	3.5	3.6	9.8	9.9	10.5	10.3
Rsym	4.4	4.5	3.6	3.7	3.1	3.6
<Intensity>	12276	8717	100706	57667	29244	23734
Scale <I> _{NEW} : <I> _{Old}	1.41		1.75		1.23	

Conclusion

The net intensities of reflections in the final data set depend on several factors:

- the amount of X-rays photons delivered by the generator in a small area
- divergence of X-ray beam
- size of crystal
- mosaicity of crystal

It is very difficult if not impossible to crystallographically separate the contributions of those individual factors. Still, by comparing the diffracted net intensities of a number of crystals of different sizes and forms, meaningful conclusions can be drawn about the performance of the X-ray equipment. In this study, we conclude that within limits of accuracy, the new GeniX ^{3D} Cu HF source is stronger than the previous model despite of running at lower power. The gain in performance for crystallographic purposes is in the same order of magnitude as the specified X-ray photon flux, up to 1.8.